

RESEARCHER RANGE

ELECTROBLOTTING APPARATUS SELECTION GUIDE

All our systems

- ✓ Have uni-directional migration.
- ✓ Cassette and electrode plate insertion in the correct direction only
- ✓ Sample loss due to user error prevented

Model	Description	Cassettes	No of gels transferred at a time	Gel size (w) x (l) (cm)	Buffer Requirement	Recommended power supply
MICROBLOT	For western blotting	One	Two	8 x 5	450 mL	Electragel-100
MIDIBLOT	For western blotting	One	Two	8 x 7	800 mL	Electracomet
TECHNOBLOT	For western, northern and Southern blotting	Two	Four	10 x 10.5	1900 mL	Electracomet

Manufactured and Marketed by **TechResource**

188, Nagdevi Street, 2nd Floor, Mumbai 400003.
www.technosource.biz, Email: ts1@vsnl.com; Tel: (022) 28937286