MIDIKIN

Standard gel size: 16.5 (w) x 16 (l) cm

- ✓ Single sided, runs one or two gels
- ✓ In-situ gel casting
- ✓ Designed for effective gel cooling
- ✓ No sample loss due to uni-directional migration
- ✓ Gel thickness: 1/ 1.5 mm
- ✓ Wells per gel: 20/ 26
- ✓ Buffer requirement: 650 mL
- ✓ 10% SDS-PAGE run in 2 hr





Features of this apparatus

Leak Proof

- Gaskets made of silicone are tough & resilient
- Leaks from plastic joints prevented as sticking is done with specially manufactured glues

Rugged

- Breakage proof PTFE combs & fibreglass notched plate
 Electrodes do not break as
- Electrodes do not bleak a they are thick & purity > 99.9%
- No corrosion due to resin filled base of male connectors

High Resolution • Sharper bands because of

- thin glass plates & deep tank
- Good cooling as entire gelsandwich is in contact with the buffer
- Low positioned electrodes
- stir buffer to cool gel

Error Minimization

- Migration of samples in wrong direction prevented due to unidirectional feature
- Wells formed straight as combs sit levelled with the
- Run up to two gels at a time
- Run up to two gels at a tim in the same system

<u>R</u>	eady-to-use system includes	05-00162 In situ casting/ Uni-directional migration	05-00163 In situ casting/ Bi-directional migration
1.	Lower anodic tank fitted with 4 mm male connector and 0.3 mm thick electrode	Y	Y
2.	Cathodic electrode unit fitted with 4 mm male connector and 0.4 mm thick electrode	Y	Y
3.	Resilient silicone tube gasket with inserts for cathodic electrode unit	Y	Y
4.	One notched and one rectangular glass plate (square bottom)	Y	Y
5.	Two 1.5 mm thick, breakage resistant, slotted spacers	Y	Y
6.	One PTFE comb (1.5 mm thick, 20 wells, est. vol. 112 µL/well)	Y	Y
7.	Set of clamps with screws	Y	Y
8.	Agarose dispenser for filling slot in spacers	Y	Y
9.	Gel-sandwich opener	Y	Y
10.	Lid for safety	Y	Y
11.	Electrical cords with male & female connectors	Y	Y
12.	Instruction manual	Y	Y

Ordering description

Midikin (05-00162): Vertical electrophoresis single sided system of gel size 16.5 (w) x 16 (l) cm, with male connectors, in situ gel casting, slotted spacers, uni-directional migration and capable of accepting two gels

Midikin-BD (05-00163): Vertical electrophoresis single sided system of gel size $16.5 (w) \times 16 (l) cm$, with male connectors, in situ gel casting, slotted spacers, bidirectional migration and capable of accepting two gels

SUGGESTION

Order along with Ready-to-use system

Gel-casting stand (05-33175), Glass plate stand (05-35165), Combs with 26 wells (05-16261/2)

Manufactured and Marketed by *TechResource*

188, Nagdevi Street, 2nd Floor, Mumbai 400003. www.techresource.in sales@techresource.in Tel: (+91) 9819095429 | (+91) 7798384170





Midikin and its spares and accessories

Item Description	Cat. No.	Unit
Midikin, ready-to-use system for uni-directional migration	05-00162	System
Midikin-BD, ready-to-use system for bi-directional migration	05-00163	System
Clamps 17 (I) cm with screws (one pair)	05-05170	Set
Screws for clamps and cathodic electrode unit (two of each)	05-06148	Packet
Resilient silicone gasket with insert for cathodic electrode unit	05-09364	Piece
Comb 1 mm thick, PTFE, 20 wells, 15 mm deep, est. vol. 75 µL/well	05-16201	Piece
Comb 1.5 mm thick, PTFE, 20 wells, 15 mm deep, est. vol. 112 µL/well	05-16202	Piece
Comb 1 mm thick, PTFE, 26 wells, 15 mm deep, est. vol. 45 µL/well	05-16261	Piece
Comb 1.5 mm thick, PTFE, 26 wells, 15 mm deep, est. vol. 66 µL/well	05-16262	Piece
Spacers, slotted, 1 mm thick, breakage resistant (one pair)	05-17171	Set
Spacers, slotted, 1.5 mm thick, breakage resistant (one pair)	05-17172	Set
Notched (3 mm) and rectangular glass plates, 20.5 (w) x 17 (l) cm (one of each)	05-28162	Set
Rectangular glass plates, 20.5 (w) x 17 (l) cm (two)	05-28172	Set
Gel-casting stand with binder clips	05-33175	Piece
Glass plate stand	05-35165	Piece
Gel-sandwich opener	61-30182	Piece
Agarose dispenser for filling slot in 1 mm spacer	62-36201	Piece
Agarose dispenser for filling slot in 1.5 mm spacer	62-36182	Piece
Electrical cords with male & female connectors design for user safety	66-40124	Pair
Gradient mixer with triple-barrel and plunging stopper. Volume per chamber 22 mL. Accommodates magnetic spin-bar of ≤ 14 mm length	67-50207	Piece
Penta-legged levelling table (23 x 31 cm) with spirit-level	69-41319	Set
Midikin, CEU. Replacement part for single sided system of gel size 16.5 (w) x 16 (l) cm, with male connector	05-01162	Piece
Midikin-BD, CEU. Replacement part for single sided system of gel size 16.5 (w) x 16 (l) cm, with male connector	05-01163	Piece
Midikin, tank. Replacement part for single sided system of gel size 16.5 (w) x 16 (l) cm, with male connector	05-02162	Piece
Midikin-BD, tank. Replacement part for single sided system of gel size 16.5 (w) x 16 (l) cm, with male connector	05-02163	Piece